

HOSPITAL



NEONATOLOGY

Complete range for the smallest patients:

- Ventilators
- Nursing care, thermotherapy and phototherapy
- Accessories

THERMOTHERAPY

LIFETHERM 2000

Warming bed

The warming bed is optionally available without cover, with closed cover or with a two-sided flip-top cover; it can be upgraded with a large selection of accessories.

- Effective heat transfer thanks to surface heating
- Precise, electronic adjustment of the target temperature
- Display of the actual temperature
- Stepless electronic height adjustment (up to 30 cm) with foot switches on both sides*
- Gentle electrical tilting of the bed surface at the touch of a button (Trendelenburg and reverse Trendelenburg positioning)
- Hinged side panels
- Two-sided flip-top or single-piece cover available
- Gel mattress
- Standard rails for mounting accessories such as catheter quiver, storage trays or suctioning accessories*
- Optimally designed for easy cleaning





Images may depict optional accessories that can be ordered separately.

Follow us:



Sales and Service:

Löwenstein Medical
Arzbacher Straße 80,
56130 Bad Ems
Deutschland
T: +49 2603 9600-0
F: +49 2603 9600-50
info@hul.de
www.hul.de

To the website:

